

Comparison of Non-radioactive Electron Capture Detector versus a Photoionization Detector for the Measurement of Diacetyl in Food

**Advanced Industrial Chemistry
5645B Paradise Blvd.
Albuquerque, NM 87114
505 890 6096
Aicmm@flash.net**

Introduction


- Objective
- Characteristics of compound
 - Why we measure it
- DBD-ECD
- Analytical set up
- Chromatograms from food stuffs
- Lessons learned/conclusions

Objectives

- Diacetyl is a compound of concern
- Test non-radioactive ECD in real-world application
 - Identify issues with detector in use
- Evaluate the option of using PID or ECD to measure diacetyl
- Both ECD and PID detectors are AIC design
- Identify issues in measuring diacetyl in foodstuffs

Diacetyl: 2,3-butanedione

- Yellow liquid with butter like odor
- Found in butter*
- Also found in beer
- Added to margarine and popcorn as flavoring



* At least that is what I had read

Why would we want to measure it?

- In beer
 - Control/monitor the beer flavor
 - Identify problems with batch
- As a flavoring
 - Margarines
 - Popcorn
- Regulatory/Health and Safety
 - NIOSH investigation of popcorn manufacturing process
 - Method 1012, Method 2557

Diacetyl, Characteristics:

- MP; -4 to -6 C
- BP; 88 C
- MW; 86 g/mole
- Density; 0.99 g/mL
- IP; 9.3 eV
- Viscinal diketone
- Kovats: DB-1 ~560, DB-Wax ~ 975

Ideal GC candidate

Two detector alternative; ECD or PID?

■ ECD:

- Advantages: highly sensitive, somewhat selective
 - Disulfides, diketones good ECD candidates
- Disadvantages: requires clean operating system (gases, no leaks), limited linear range, **radioactive source**

■ PID:


- Advantages: highly sensitive, somewhat selective; only one gas required
- Disadvantages: limited linear range, not as selective as ECD, lamp window may need cleaning.

One better than the other?


DBD-ECD eliminates biggest disadvantage of ECD detector: replace radioactive source

- DBD = Dielectric Barrier Discharge plasma
 - AC discharge across a dielectric barrier
 - Non-thermal discharge
 - Low electrode wear/large electrode surface
 - Ability to operate without getters/purging
- Simple design
 - Non-radioactive, windowless
 - Simple, robust power supply
 - Conventional electrometers


DBD-ECD Schematic:



Detector Picture



Installed on Varian 3400



Helium plasma color

Analytical set up: extraction

- Headspace Extraction (beer, butter, margarine)
 - Sample material into 20 mL VOA vial
 - 3 grams butter or margarine, 15 mL beer
 - Heated for 10 min at 50C in GC oven
 - 4 mL room air into vial, 3 mL vial air out
 - Aliquot transferred to GC using warm syringe
- Popcorn
 - Un-popped and microwave popped
 - From the bag, syringe into bag to withdraw sample
- Standards prepared in warmed static dilution bulb

Analytical set up: instrumental

- HP 5890 A
 - 6-Port sample valve, 0.1 mL loop (65C)
 - Dual column configuration; Valco T inside oven
 - 30 meter RTX-VMS, 0.25 id X 1.4 film to ECD
 - 30 meter RT-Stabiliwax, 0.25 id X 0.5 film to PID
 - Flow rate: 6 mL/min (split to both columns)
 - Oven program: 55C/6 min/14C/min/95 C/1 min
- Dual column configuration
 - Allows for confirmation (unsure of early interferences)
 - Real time, direct comparison of detectors (same sample to both)
 - Allows evaluation of ECD selectivity

Instrumental set up; Dual Detectors


■ DBD-ECD

- Non-radioactive source (dielectric barrier discharge)
- Uses helium reaction gas (five 9's), hydrogen dopant
- Copper plumbing, no extra purification, typical GC flow controllers
- Stacked electrode configuration
- Uses HP 5890 ECD electrometer (N2 mode)


■ In-house PID Detector

- 10.2 eV lamp
- Stacked electrode configuration
- Use HP 5890 FID electrometer
- Hydrogen sweep gas (10 mL/min)

Standard chromatogram both channels; 800 pg/column




Butter (3 gram extraction)




DA not detected (PID); very low on ECD*

Margarine; both channels



Much higher levels than butter

Margarine; ECD channel




No significant interferences

Comparison of two margarines and agreement between ECD and PID


	DA PID	DA ECD	PID/ECD	ng/g (ECD)
M1	5.89	7.07	0.83	58.88
M1 xt	9.72	10.29	0.94	85.78
M2	13.97	13.32	1.05	111.02
M2 xt	18.37	16.42	1.12	136.81

Chromatogram from freshly opened stout




Even ECD subject to some interferences


Chromatogram from flat beer co-injection, both channels: ECD clear of interferences



Chromatogram from flat beer co-injection: PID channel subject to interference



Chromatogram from freshly popped popcorn (PID)



Significantly less on second analysis (sample cooled)

Lessons learned or,

“If you want to do Diacetyl in food...”

- **Headspace optimization advisable**
 - Longer incubation yields more sensitivity
 - Carbonation not an issue under these conditions
- **Need a better way to sample popcorn**
 - Sample portion of a bag
 - Should be microwave safe
 - Sample very wet directly after popping
 - Sampling immediately after popping
- **Dual column configuration advantageous**
 - Especially if using PID's or FID's
 - Subject to interferences (especially ethanol)

Conclusion, Diacetyl in food stuffs

- Clearly evident in margarine, beer, popped popcorn
- Levels much lower than expected when starting out
 - Could be a function of headspace conditions
 - Lower levels in the flat beer than the fresh beer
- Very low concentration in butter (??)
 - Much higher levels in margarine than butter
- Difference in levels between two margarines
 - But is it significant?

Conclusion, ECD vs. PID

- **Both ECD and PID can be used to measure Diacetyl**
 - **Depends on the levels you need to measure**
- **ECD**
 - **Much, much better detection limit**
 - **More selectivity**
- **PID**
 - **Sensitive enough for some samples (margarines)**
 - **Subject to interferences (i.e. ethanol)**
 - **Able to detect other constituents (see beer and popcorn)**
- **Good agreement between two detectors when free of interferences**

Conclusion, DBD-ECD

- **Dielectric Barrier Discharge can be used as a source for ECD**
 - **Stable plasma source, using five 9's helium**
 - **Does not require getters or specialized plumbing**
 - **Uses conventional ECD electrometers**
- **DBD-ECD very sensitive to vicinal diketones**
- **DBD-ECD selective**
 - **Absence of significant ethanol peak in beer analysis**
- **DBD-ECD non-rad source eliminates ownership issues**
 - **No licensing, swiping, custody issues**
 - **Can be easily dis-assembled and cleaned**