

Saving Helium on the 5890, 6890 and 7890 GC

Matthew Monagle

Advanced Industrial Chemistry LLC

4804 2nd Street

Albuquerque NM, 87107

505 890 6096

aicmm@flash.net

Why look at helium switching?

- Finite resource:
 - Save it for essential applications like Helium Ionization Detectors (HID's)
- Supply concerns
 - Sometimes difficult to acquire;
- COST, COST, COST!
 - Over \$400/bottle

Conserve the expensive gas for necessities

Alternatives to carrier switching?

- Switch permanently; move to hydrogen
 - Excellent chromatography
 - Concerns about flammability
- Move to some other carrier;
 - Chromatography or run time may suffer
 - Detector sensitivity (e.g. TCD) may suffer

Where possible, H₂ is an excellent alternative

Why stay with helium?

- Good carrier;
 - Excellent chromatography
 - Safe
- Some instruments (i.e. HID's) require it
- Some methods require it
 - e.g. Numerous ASTM methods
 - Detector sensitivity (e.g. TCD) may require it
 - Non-reducing carrier (HT SimDis)


When you have to stay with helium...

Switching is way to stay with helium


- Move to helium to run method specified helium analysis;
 - GC control file switches to helium
 - Wait long enough to clear out idle gas (blank run)
- Move back to a different carrier for idle times
- Only consuming helium during analytical operations

Simple way to have the best of both


Implementation: Switching valve is simple installation



Hardware configuration



5890



7890

Controlled by GC software

Procedure:

- Set up sequence with awake method
 - Switches from idle gas (H₂, N₂, etc...) to helium carrier
- Runs a run long enough to effect carrier gas changeover
 - Time established empirically
- Runs samples (some number)
 - Regular analytical method
- Runs a run to switch back to alternative carrier (sleep method)

Sequence Example

Sequence Table: D3606

Currently Running

Line:

Method:

Location:

Inj:

Sample Info

--

Line	Location	Sample Name	Method Name	Inj/Location	Sample Type	Ca
1			AWAKE	1	Sample	
2	Vial 1	QC	BENZENE	1	Sample	
3			SLEEP	1	Sample	

Method screen shot (5890) Analytical/Awake

File RunControl Instrument Method Sequence View Abort Help

Valves/Relays: Instrument 1

Initial Setpoints

5890 Valves On: Valve 1 Valve 2 Valve 3 Valve 4

19405 Relays On: Relay 1 Relay 2 Relay 3 Relay 4

19405 Relays On: Relay 5 Relay 6 Relay 7 Relay 8


Valve/Relay Time Table

Time	Name	State	Comment
5.20	Valve 1	<input checked="" type="radio"/> On	
	Valve 2	<input type="radio"/> Off	
	Valve 3		
5.20	Valve 2	ON	

Buttons: Add, Replace, Copy, Paste, Cut, OK, Cancel, Help


Switching valve on

Method screen shot (5890) Sleep



Switching valve off

3606 Chromatogram (includes switching)



Implemented on:

- 5890:
 - 2 instruments, manual pneumatics, helium/nitrogen, 1 TCD, 1 FID, 1 packed column, 1 capillary column
- 6890:
 - 1 Instrument, EPC pneumatics, helium/hydrogen, FID, capillary column
- 7890
 - 1 Instrument, EPC pneumatics, helium/hydrogen, FID, capillary column

Where it is not appropriate:

- Sensitive detectors
 - Cannot recover fast enough to trace amounts of idle gas (e.g. HID's)
 - Not “bubble tight”; might see air contamination
 - Mass Spectrometer: hydrogen switch over can be very dirty
- Applications that might measure sleep gas
 - Some TCD applications not appropriate

Where it is not appropriate (cont.):

- Method used all day, day in and day out
 - Not worth switching; hard on the switching valve
 - Not really cost effective
- May not be appropriate for certain instrument combinations
 - In my opinion: Switching to hydrogen idle with a TCD (safety issue) with a 7890

When you do use switching, be
aware . . .

Identified issues:

- 7890 EPCs
 - Cannot re-configure; use hydrogen
 - Instruments converted are FID's so flammability not an issue
- 5890 manual pneumatics
 - No EPC - use nitrogen – cheap, safe
- Default to valve off (saves wear and tear)
 - Must return to sleep

Cost Savings:

- Working with customer to ascertain helium expenses for 2012 vs. 2013
- Also looking at H₂ consumption for 2012 vs. 2013
- Complicated by move to H₂, N₂ generators

When you have to stay with helium carrier ...

Conclusion:

- Carrier switching offers a simple inexpensive means of reducing helium consumption in the laboratory
- Analytical methods stay with helium carrier
- Implementation is simple and inexpensive
- Not applicable to all instruments, all applications